







Search



Home Products Supports Company Contact Us

Products / Physics Lab Instruments / Polarization / Brewster's Angle Apparatus

Spectroscopy

Interferometry

Polarization

Brewster's Angle Apparatus

Malus Law Apparatus Faraday Effect Apparatus Kerr Effect Apparatus Pockel Effect Variable Angle Laser Ellipsometer

Fiber Optics

Lasers & Opto-Electronics

Diffraction

Optics

Electricity & Magnetism

Acoustics

Holography

Sources & Detectors

Brewster's Angle Apparatus

Model No: HO-ED-P-01


A beam of light incident on a dielectric transparent material can be resolved into parallel (P) and orthogonal (S) components. These components have different reflection coefficients and Brewster discovered that at a particular angle of incidence $\theta \mathtt{B}$ (called Brewster angle), reflection co-efficient of P-component goes zero. At this angle direction of reflected and transmitted beam are orthogonal to each other.




Experiment Examples

- ▶ To measure and plot the graph-reflectivity versus angle of incidence.
- ▶ To find the Brewster's angle (also known as the polarization angle) of glass plate and determine it's refractive index.

This equipment uses goniometer and a pinhole photo detector to determine the Brewster's angle and to study polarization of reflected light. Intensities of reflected light polarized in the plane of incidence and perpendicular to the plane are obtained as a function of angle of incidence. The results should be consistent with Fresnel's laws of reflection. All components are made out of anodized aluminium and stainless steel for corrosion free life.





Scope of Supply

Goniometer with Detector Mount

Model No: ED-P-01-GDM

Kinematic Laser Mount

Accessories

Resolution of prism table 10 arc min Resolution of stage 1 arc min Quantity



Optical Rail

Polarizer Rotator with Mount	+

Glass Slide			+

Diode Laser with Power supply (Red)	+

Detector Output Measurement Unit	+



Home | Products | Supports | Company | Contact Us Privacy Policy | Terms & Conditions © 2014 Holmarc Opto-Mechatronics (P) Ltd. All rights reserved.





